

Why choose the X-COLLAR?

Until now, the marketplace has not offered a cervical collar with sufficient impediments to adequately immobilize the head and cervical regions of a trauma patient. In fact, if one were to put on a commonly used cervical collar one would find insufficient immobilization of the head and neck regions which can easily flex and extend, as well as, turn from side to side. According to our research, most cervical collars allow up to 90 degrees of free movement from side to side. This lack of support and immobilization could worsen existing injuries.

Another troublesome deficiency of today's cervical collars is that they were not designed to be adjusted to the patient after application. Currently, the collar size must be approximated and adjusted prior to application. If the adjustment is incorrect the collar must be removed, re-adjusted and re-installed. The emergency medical field is in need of a cervical collar that will better immobilize the patient's head, neck, and thoracic regions while being fully adjustable.

The X-COLLAR is the answer. It is the only custom fitting cervical immobilization device that meets the demands of today's emergency medical field:

- Enhanced immobilization
- Adjusts after application
- Immobilizes patients in position of comfort
- Custom fits a broader range of patients



Martin Rizo Patron

Phone
805-637-6492
Fax
805-566-1403
mrizopa@hotmail.com
www.xcollar.com

X COLLAR

A NEW STANDARD OF CARE



X COLLAR

RIZO PATRON ENTERPRISES

designed by emergency care professionals

RIZO PATRON ENTERPRISES

X-COLLAR

The X-COLLAR is designed by emergency care professionals to address the common deficiencies and disappointments of other cervical immobilization devices by featuring:

- Enhanced immobilization
- Adjusts after application
- Immobilizes in position of comfort
- Custom fits a broader range of patients


ENHANCED IMMOBILIZATION

Protecting the patient from movement that may cause further injury is the X-COLLAR'S primary objective.



Inadequate immobilization is one of the biggest downfalls of other collars which can allow the head up to 90 degrees of side to side movement.

The X-COLLAR is designed to protect and support the cervical region by securing it to the patient's torso for maximum immobilization. This is largely due to the collar's enhanced back and chest supports.



ADJUSTABILITY

Post application bi-lateral and vertical adjustments allow the X-COLLAR to better fit the patient's body type while reducing manipulation of the spine.

Products currently used in the field force the patient's neck into the collar's shape creating additional discomfort and manipulation of the damaged spine.

POSITION OF COMFORT

The exclusive dual-sided adjustable straps allow the X-COLLAR to immobilize the patient's head in a variety of positions.



RANGE

The X-COLLAR is so flexible that it may be fitted from a small child to a large adult.

28 Vertical adjustments (3.25in.)
56 bilateral adjustments (5.5in.)

Color-coded numbers allow precise adjustments



The pushbutton locks the chest support in place.



Color-coded velcro-straps make application quick and easy.

